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# CHONDROreparant®

**HYALREPAIR®-10 Chondroreparant®**  
(how supplied – a glass vial, 5.0 ml)

Gel based on modified hyaluronic acid hydrous sterile

Total content of salts of hyaluronic acid – 0.8%

## Usage instruction

### Ingredients

Zinc hyaluronate and copolymer of sodium hyaluronate with sodium ascorbyl phosphate – 8 mg/ml, L-cysteine – 1 mg/ml, L-glutathione – 3 mg/ml, sodium chloride, phosphate buffer, water for injection.

### Description

HYALREPAIR®-10 Chondroreparant® is a transparent, colorless, homogeneous, gel-like substance, supplied in a glass vial.

Contains a unique form of high-purity modified hyaluronic acid obtained by bacterial fermentation with addition of hyaluronic acid copolymer with sodium ascorbyl phosphate, L-cysteine and L-glutathione.

### Mechanism of action

The state of fasciae and tendons, viscoelastic characteristics, defensive and trophic functions depend largely on concentration of specific proteoglycan – hyaluronic acid, vitamin C and amino acids, which provide dynamic balance between biosynthesis and catabolism of collagen and elastin.

Joint immobilization is accompanied by morpho-functional changes in locomotor system tissues.

Excessive physical exercises, traumas, hypoxia, microcirculation and metabolic tissue disorders result in excessive adnation of a tendon with surrounding tissues, degenerative disorders with increased formation of covalent bonds between collagen molecules with formation of tendogenic and myogenic contracture and their transformation to fibrous (cicatrical) tissue in serofibrous and chronic fibrous tendinitis. Elevated levels of homocysteine and dimethylarginine in matrix of fascial tissue form aseptic inflammation, which results in ruptures and other dama-

ges of tendon fascicles, adipose tissue, nerves, lymphatic and blood vessels.

Oxidation-reduction system malfunction is one of the reasons of microcirculation abnormality. Increased production of reactive oxygen results in ischemia (hypoxia) most noticeable in myofascial tissues, containing furcated capillary network.

HYALREPAIR®-10 Chondroreparant® is a bioactive composition of modified hyaluronic acid, ascorbic acid, L-glutathione and L-cysteine which participate in oxidation-reduction processes and in regulation of collagen and elastin initial synthesis. It makes supplies of the complex and allows targeted delivery of small doses of active components to damaged connective tissue formations. Being antioxidants, glutathione, cysteine and ascorbic acid inhibit pathologic changes caused by toxic metabolites and xenobiotics in all types of connective tissue and cells.

HYALREPAIR®-10 Chondroreparant® action is based on biologic stimulation of myotendinous normal metabolic processes which improves physiologic and rheologic status of trophic, growth regulation and cell differentiation systems.

Structural hyaluronic acid modification to HYALREPAIR®-10 Chondroreparant® provides prolonged reparative action on tissues damaged due to traumas or degenerative-dystrophic processes. The maintenance of optimal state of intercellular matrix of cartilage and tendinous-ligamentous apparatus also depends on the glycosaminoglycan synthesis and degradation balance.

One of the components, vitamin C, participates in collagen triple-helix stabilization and in collagen molecules covalent bonds building for collagen fibrils formation. Ascorbic acid decreases the excess urates formation in connective tissue dystrophic calcification.

Intracellular L-glutathione inhibits production of inflammation mediators. L-cysteine, being a donor of sulfhydryl groups, takes part directly in recapping of free radicals whereby discontinuing free-radical oxidation reactions and defending cellular information centre.

Zinc mediates synthesis and degradation of various neurotransmitters reducing pain sensation. It accelerates wound healing due to phagocyte formation and macrophages activity increase, contributing to necrosis elimination, activates collagen necrosis elimination, activates collagen synthesis mediating enzymes and attracts fibroblasts to the damage area. It is very important for antibacterial immunity. It fulfills hyaluronic acid deficiency and activates matrix metalloproteinases (MMPs), which are able to disintegrate almost all components of extracellular matrix of connective tissues.

### Indications for use

HYALREPAIR®-10 Chondroreparant® is used in sports medicine, rheumatology, traumatology, orthopedics and maxillofacial surgery.

Periarticular, intramuscular and perineural administration is allowed.

The product is used:

1. As part of complex therapy for:

- periarthrit, tendinitis, tenosynovitis stenosans (de Quervain disease), bursitis, enthesitis;
- radicular syndrome, spinal stenosis, fibrotic changes of muscles (proliferative myositis), ligamenti, fascial restrictions (proliferative fasciitis), periarticular algetic fibrotic changes, trigger points of repetitive stress injury;
- subfascial hypertensive syndrome, tendinopathy, enthesopathy;

- spondyloarthrosis and facet syndrome, diffuse idiopathic skeletal hyperostosis (Forestier disease), ankylosing spondylitis (Bekhterev's disease);
- degenerative-dystrophic changes of ligamentous apparatus, connective tissue formations and cartilages;
- inactive periarthrit.

2. During the rehabilitation period after injuries and surgical interventions.

### Contraindications

1. Hypersensitivity to the components.
2. Anamnesis of autoimmune diseases or current autoimmune therapy.
3. Bleeding sickness (endogenic or secondary to anticoagulants use).
4. Infectious (septic) inflammatory process of joints or periarticular tissues, general infectious disease.
5. Signs of active skin disease or skin infection near the injection site.
6. Pregnancy and lactation, children.

### Precautions

Do not use disinfectants containing quaternary ammonium salts. Before the injection it is necessary to warm the gel up to the room temperature.

Patients should be warned about restriction of alcohol drinking during all treatment course by HYALREPAIR®-10 Chondroreparant®. Ethanol increases the amount of free radicals and disturbs intracellular metabolism.

After the injection of HYALREPAIR®-10 Chondroreparant® the patient is recommended avoid excessive physical activity for several days.

Do not use HYALREPAIR®-10 Chondroreparant® in the same syringe with the anesthetic and other pharmacological drugs.

There have been no investigations done on the results of HYALREPAIR®-10 Chondroreparant injections to women during pregnancy and lactation and children.

### Administration method

Dosage and deepness of injection of HYALREPAIR®-10 Chondroreparant® to projection of damage area of ligamentous apparatus and tendons, as well as selection of injection needle in each case is within the discretion of a professional.

Depending on anatomic-physiologic characteristics and soft tissues state multiple injection of 0.2-0.5 ml fortnightly during 6-8 weeks is recommended. The treatment course should be repeated after 6-12 months. Concomitant injections in several zones are possible.

**Attention!** To achieve an effective result of the treatment it is vitally important to perform the injection correctly.

### Not to be used for intra-articular injection!

The injection of HYALREPAIR®-10 Chondroreparant® may be performed only by professionals having a special permission according to the local regulations. "Laboratory THOSCANE" LLC does not take responsibility for the consequences of the incorrect use of HYALREPAIR®-10 Chondroreparant®.

### Side effects

Periarticular, intramuscular and perineural injection of HYALREPAIR®-10 Chondroreparant® is well tolerated. In rare cases secondary local events are possible, such as pain, feeling of warmth, effusion and extremely rarely allergic reactions.

### Interaction with other medical devices

No data concerning interaction of HYALREPAIR®-10 Chondroreparant® with other formulations for extra-articular injection.

### Disposal

Used gel should be disposed in compliance with local regulations.

**The vial, the syringe, the needle and the remaining unused gel should be disposed of immediately after the procedure.**

### How supplied

5.0 ml of HYALREPAIR®-10 Chondroreparant® is supplied in a dark glass vial. 1 needle 27G ½" (0.4×13 mm) is included into a set. The content is sterile. A glass vial with 2 adhesive labels for physician attached is packaged with instruction for use and needle in a cardboard box. Apply adhesive label for physician to a patient's case history to ensure future gel identification.

### Storage terms and shelf life

Store in a place protected from sunlight at the temperature between +5 and +25 °C. Do not freeze. Do not use after expiration of indicated shelf life.

**DO NOT USE THE PRODUCT IF THE PACKAGE IS DAMAGED!**

Complaints and claims on quality of medical device are accepted by an authorized representative of the manufacturer.

### Manufacturer

"Laboratory THOSCANE" LLC, Russia

### Symbols on the package

Caution!  
Consult instructions for use

Do not reuse

Temperature limitation

Keep away from sunlight

Sterilized using steam or dry heat

Sterilized using irradiation

Do not resterilize

Do not use if package is damaged

Use by date

Batch code

Manufacturer