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
## HYALREPAIR®-02 Chondroreparant®

(how supplied – a glass syringe, 2.0 ml)

Gel based on modified hyaluronic acid  
hydrous sterile

Total content of sodium hyaluronate – 1.5%

### Usage instruction

#### Ingredients

Sodium hyaluronate and copolymer of sodium hyaluronate with magnesium ascorbyl phosphate – 15 mg/ml, L-proline – 3 mg/ml, L-lysine hydrochloride – 3 mg/ml, glycine – 3 mg/ml, sodium chloride, phosphate buffer, water for injection.

#### Description

HYALREPAIR®-02 Chondroreparant® is a transparent, colorless, homogeneous, gel-like substance, supplied in a glass syringe.

Contains a unique form of high-purity modified hyaluronic acid obtained by bacterial fermentation with addition of hyaluronic acid copolymer with magnesium ascorbyl phosphate, L-proline, L-lysine and glycine.

#### Mechanism of action

The state of metabolism in joints and viscoelastic and protective properties of the synovial fluid are considerably influenced by the concentration of hyaluronic acid in it, a specific proteoglycan found in cartilage, synovial membrane and other connective tissue structures. The action of HYALREPAIR®-02 Chondroreparant® is based on the biological stimulation of natural metabolic processes, which results in the improvement of physiological and rheological status of the synovial fluid.

HYALREPAIR®-02 Chondroreparant® reduces the rate of destruction of articular cartilage and stimulates the restoration of the synovial fluid homeostasis.

As a biologically active combination of modified hyaluronic acid with ascorbic acid and amino acids HYALREPAIR®-02 Chondroreparant® takes part in the collagen synthesis and allows for targeted delivery of the active substances directly to the damaged connective tissue structures.

Structural modification of hyaluronic acid in HYALREPAIR®-02 Chondroreparant® determines its prolonged reconstructive action on the cartilages, ligaments or connective tissue structures of the joint damaged as a result of injuries or degenerative-dystrophic processes.

Maintenance of the optimal condition of the intercellular matrix of cartilage and tendinous-ligamentous apparatus also depends on the balance of synthesis and degradation of this glycosaminoglycan. Following the intra-articular injection HYALREPAIR®-02 Chondroreparant® releases amino acids including essential acid lysine. Together with vitamin C they create physiologically favorable environment for the enhancement of the cells' metabolic activity

and stimulate the synthesis of the main components of the intercellular matrix of cartilage and connective tissue structures.

HYALREPAIR®-02 Chondroreparant® restores the effective concentration of hyaluronic acid in the matrix which results in stimulation of fibroblasts and chondrocytes. Intermediate metabolites of the hyaluronic acid degradation produce angiogenic, immunostimulating, and anti-inflammatory effects.

Vitamin C is used in metabolism of glycine, proline and lysine. Proline hydroxylation is necessary for stabilization of the collagen triple helix and lysine hydroxylation is very important for further covalent bonding between molecules of collagen in the process of collagen fibrils formation. This process is essential for stabilization of the intercellular matrix of all types of connective tissue.

#### Indications for use

HYALREPAIR®-02 Chondroreparant® is used in sports medicine, rheumatology, traumatology, orthopedics and maxillofacial surgery. It may be administered into any joints.

Indications:

1. As a part of the combined therapy:
  - treatment of I-III stage osteoarthritis of knee, coxofemoral and other synovial joints of degenerative-dystrophic or traumatic origin;
  - necessity to improve viscoelastic and protective properties of the synovial fluid;
  - degenerative-dystrophic conditions of the ligamentous apparatus, connective tissue structures and cartilages;
  - arthritis in remission.
2. For prevention of degenerative-dystrophic diseases of the musculoskeletal system.
3. During the rehabilitation period after injuries and surgical interventions.

#### Contraindications

1. Hypersensitivity to the components.
2. Anamnesis of autoimmune diseases or current autoimmune therapy.
3. Bleeding sickness (endogenous or secondary to anticoagulants use).
4. Severe inflammation or infection of the joint.
5. Signs of active skin disease or skin infection near the injection site.
6. Pregnancy and lactation, children.

#### Precautions

It is necessary to take general precautions for intra-articular injections.

HYALREPAIR®-02 Chondroreparant® must not be injected in the presence of infected wounds, scratches in the joint area, evident inflammation, significant effusion and venous or lymphatic stasis on the side of the damaged joint.

Do not use disinfectants containing quaternary ammonium salts. Before the injection it is necessary to warm the gel up to the room temperature.

HYALREPAIR®-02 Chondroreparant® must be injected accurately into the joint cavity.

After the injection of HYALREPAIR®-02 Chondroreparant® the patient should observe a relaxed regimen and avoid excessive exertion for several days.

Do not inject anesthetics or other pharmacological preparations into the affected joint cavity during the course of HYALREPAIR®-02 Chondroreparant® injections.

There have been no investigations done on the results of HYALREPAIR®-02 Chondroreparant® injections to women during pregnancy and lactation and children.

#### Administration method

The amount of the injected gel and the type of a needle are determined by a physician on a per-patient basis depending on the anatomic and physiological characteristics and the condition of the joint.

In case of effusion removal of the excessive synovial fluid from the joint before the injection of the gel is recommended.

Depending on the size of the joint, it is necessary to inject into the joint up to 2 ml of HYALREPAIR®-02 Chondroreparant® every 14 days for 6-8 weeks. Refresher course is conducted in 6-12 months. It is possible to treat several joints simultaneously.

Treatment with HYALREPAIR®-02 Chondroreparant® produces the effect only on the joint into which the gel is injected, HYALREPAIR®-02 Chondroreparant® does not produce a systemic effect.

#### Attention!

To achieve an effective result of the treatment it is vitally important to perform the injection correctly. The injection of HYALREPAIR®-02 Chondroreparant® may be performed only by professionals having a special permission according to the local regulations.

“Laboratory THOSCANE” LLC does not take responsibility for the consequences of the incorrect use of HYALREPAIR®-02 Chondroreparant®.

#### Side effects

Intra-articular injection of HYALREPAIR®-02 Chondroreparant® is well tolerated. In rare cases secondary local events are possible, such as pain, feeling of warmth, effusion and extremely rarely allergic reactions.

#### Interaction with other medical devices

Introduction of HYALREPAIR®-02 Chondroreparant® in combination with other medical devices has not been researched.

#### Disposal

Used gel should be disposed in compliance with local regulations.

**The syringe, the needle and the remaining unused gel should be disposed of immediately after the procedure.**

#### How supplied

2.0 ml of HYALREPAIR®-02 Chondroreparant® is supplied in a glass syringe with Luer-Lock tip. The content is sterile. Syringe with gel, 2 adhesive labels for physician attached to a syringe are enclosed in a sealed blister. One blister is packaged with instruction for use in a cardboard box. Apply adhesive label for physician to a patient's case history to ensure future gel identification.

#### Storage terms and shelf life

Store in a place protected from sunlight at the temperature between +5 and +15 °C. Do not freeze. Do not use after expiration of indicated shelf life.

#### DO NOT USE THE PRODUCT IF THE PACKAGE IS DAMAGED!

Complaints and claims on quality of medical device are accepted by an authorized representative of the manufacturer.

#### Manufacturer

“Laboratory THOSCANE” LLC, Russia

#### Symbols on the package

Caution!  
Consult instructions for use

Do not reuse

Temperature limitation

Keep away from sunlight

Sterilized using steam or dry heat

Do not resterilize

Do not use if package is damaged

Use by date

Batch code

Manufacturer